CHEMI Ca L REACTIONS AND EQUATIONS

In this experiment, you will observe changes that occur during chemical reactions. You will also balance the corresponding chemical equations and identify the reaction type. You will also get familiarized to chemical nomenclature while you encounter chemical names and formulas as you work.

In chemical reactions, a set of substances is transformed into another set of substances with different chemical identity. For example, silver tarnishes because it reacts with hydrogen sulfide in air. The chemical reaction taking place is indicated by a change in color. The shiny silver metal turns dull dark-gray because a new is substance, silver sulfide, is formed on the surface. In addition to a change in color, other evidences of a chemical reaction occurring are formation of gas (bubbles or fizzing), formation of a precipitate (insoluble solid) and heat given off or absorbed (change in temperature).

During chemical reactions, atoms of the reactants rearrange into new combinations to form products. By the law of conservation of matter, the total number of atoms of each element in the reactants is equal to the total number of atoms of each element in the products. A chemical reaction is concisely represented by a balanced chemical equation.

In a chemical equation, the reactants are shown on the left side and the products on the right, separated by an arrow. When silver tarnishes, the reactants are solid Ag and gaseous H₂S, and the products are solid Ag₂S and gaseous H₂. To balance the equation, coefficients must be used so that the numbers of atoms of each element are the same on both sides of the arrow:

Unbalanced Equation:	$\mathrm{Ag}(s) + \mathrm{H}_2 \mathrm{S}(g) \rightarrow \mathrm{Ag}_2 \mathrm{S}(s) + \mathrm{H}_2(g)$
Balanced Equation:	$2 \operatorname{Ag}(s) + 1 \operatorname{H}_2 S(g) \rightarrow 1 \operatorname{Ag}_2 S(s) + 1 \operatorname{H}_2(g)$

Most chemical reactions can be classified into the following types:

1. Combination or Synthesis – **Two or more** simpler substances (compounds or elements) form **one** complex substance

2. Decomposition – **One** complex substance is breaks up into **two or more** simpler substances (compounds or elements).

General Format:	С	\rightarrow	А	+	В
Example:	$2 \operatorname{NaN}_{3}(s)$	$\rightarrow 2$	Na (s)	+	$3 N_2(g)$

3. Single displacement – Reaction between **an element and a compound** where one element takes the place of another in the compound.

General Format: A + BY \rightarrow B + AY Example: 1 Ca (s) + 2 HCl (aq) \rightarrow 1 CaCl₂ (aq) + 1 H₂ (g)

4. Double displacement – Reaction between two compounds in which ions switch partners.

 Combustion – Substance reacts with oxygen to form an oxide product. (Note: Reacting a hydrocarbon (e.g. CH₄) with O₂ always produces CO₂ and H₂O, both of which are technically oxide products)

*The reaction of a metal with oxygen forming an oxide can also be classified as combination/synthesis.

PROCEDURE

For each of the procedure below, record your observations. Remember that you are looking for the following evidences that a reaction has occurred:

- New substance formed (change in color)
- Formation of gas (bubbles or fizzing)
- Formation of a precipitate (insoluble solid)
- Heat given off or absorbed (change in temperature)
- Light emitted

Consult the chemical equation in the worksheet to confirm your observation is consistent with the expected products. Balance the chemical equation for the reaction that occurs. Identify the reaction type.

For Parts I, III, and V, use small quantities and note that the quantities need not be exact. For solids, use the amount that fills the tip of a spatula. For liquids, use ~ 3 mL. (Tip: Use a graduated cylinder to measure 3 mL water. Transfer to a test tube to see and estimate about how much liquid to use each time the procedure calls for 3 mL).

IMPORTANT: Dispose wastes as soon as you finish a reaction or part so you can keep track of the test tube contents and place them in their designated containers.

I. Zinc and Copper(II) Sulfate

- 1. Place ~3 mL 1 M solution of copper (II) sulfate (CuSO₄), into each of two test tubes.
- 2. Obtain a small piece of zinc metal.
- 3. Observe and record the appearance of both CuSO₄ solution and zinc metal.
- 4. Add the zinc to the CuSO₄ solution in one test tube. The CuSO₄ in the other test tube serves as reference for initial color.
- 5. Observe after 15 minutes. Observe again after another 15 minutes (30 minutes total).

II. Magnesium and Oxygen

- 1. Obtain a small strip of magnesium ribbon and make initial observations.
- 2. Hold the end of the magnesium ribbon with tongs. Ignite it using a Bunsen burner flame.

CAUTION: Remove the magnesium from the flame as soon as it ignites. Do not look directly at the flame because the bright light can damage your eyes.

3. Record the changes that occur.

III. Metals and Hydrochloric Acid

1. Place \sim 3 mL 1 M solution hydrochloric acid (HCl) in each of three separate test tubes.

CAUTION: HCl is corrosive. Clean up spills immediately. If it spills on the skin, put affected area under running water for at least ten minutes.

- 2. Obtain a piece of copper and zinc. Record the appearance of the metals.
- 3. Add one metal to the acid in each of the test tubes.
- 4. Record any changes that occur.
- 5. Obtain a piece of magnesium and a wooden splint. Record the appearance of the metal.
- 6. Ignite the wooden splint, then immediately place the magnesium into the third test tube with HCl. Wait a few seconds, then insert the burning wooden splint into the top of the test tube (do not place it in the liquid).

CAUTION: Dispose of the wooden splints in the designated waste container. Do not place it in the trashcan.

IV. Reactions of Ionic Compounds

The procedure that follows is for mixing 0.1 M solutions of:

- A. calcium chloride (CaCl₂) and sodium phosphate (Na₃PO₄)
- B. barium chloride (BaCl₂) and sodium sulfate (Na₂SO₄)
- C. iron (III) chloride (FeCl₃) and potassium thiocyanate (KSCN)
- 1. Record initial observations for each solution you are mixing.
- 2. Obtain three separate test tubes, one for each pair of solutions (A, B and C) above.
- 3. To each test tube, add 20 drops of both solutions. Record any changes that occur. Watch for changes in color and formation of gas or solid.

V. Sodium Carbonate and Hydrochloric Acid

1. Place \sim 3 mL of 1 M HCl in a test tube.

CAUTION: HCl is corrosive. Clean up spills immediately. If it spills on the skin, put affected area under running water for at least ten minutes.

- 2. Add 0.5 g of solid sodium carbonate (Na₂CO₃) to the test tube. Record your observations.
- 3. Light a wood splint and insert the flame in the neck of the test tube. Note what happens to the flame. Explain your observation.

CAUTION: Dispose of the wooden splints in the designated waste container. Do not place it in the trashcan.

CLEAN-UP

- Dispose of wastes in designated containers in the front hood.
- Place wooden splints in the designated waste container. Do not place it in the trashcan.
- Wash all glassware used. Return materials where they belong. Shake off excess water from washed test tubes and return to the test tube racks on your work station.
- Make sure the Bunsen burner and gas valve are off, disconnect the tubing from the gas valve.

Date:	
-------	--

Name: ______
Partner's Name: ______

CHEMICAL REACTIONS AND EQUATIONS

DATA AND OBSERVATIONS

When balancing the equations, write down all coefficients, including 1. Use the appropriate terms (see pg 2) for describing the evidence of a chemical reaction, include colours where appropriate.

I. Zinc and Copper(II) Sulfate

	Observations			
Time elapsed	Solid	Solution		
Initial				
After 15 min				
After 30 min				
	Metal formed:	Solution formed:		

Balance: $__Zn(s) + __CuSO_4(aq) \rightarrow __Cu(s) + __ZnSO_4(aq)$

II. Magnesium and Oxygen

Observations:

Initial _____

During Reaction	
0	

Balance: $Mg(s) + O_2(g) \rightarrow MgO(s)$

Reaction type: _____

III. Metals and Hydrochloric Acid

	Observations		
	Cu	Zn	Mg
Initial			
Reaction with HCl			
Balance: Cu(s) +	$\underline{\qquad} \operatorname{HCl}(\operatorname{aq}) \to \underline{\qquad}$	CuCl ₂ (aq) + H ₂ (g)
Zn(s) +	$- \underline{\qquad} \operatorname{HCl}(\operatorname{aq}) \rightarrow \underline{\qquad}$	$\underline{ZnCl_2(aq)} + \underline{H_2}$	(g)
Mg(s) +	$- \underline{\qquad} \operatorname{HCl}(\operatorname{aq}) \rightarrow \underline{\qquad}$	$_{\rm MgCl_2(aq)} + _{\rm H_2}$	e(g)
Reaction type:			
What caused the noise	when the splint was in	nserted into the test tub	e?

IV. Reactions of Ionic Compounds

	Observations			
	A. CaCl ₂ and Na ₃ PO ₄	B. $BaCl_2$ and Na_2SO_4	C. FeCl ₃ and KSCN	
Initial				
Reaction				

Balance:

A. _____CaCl₂(aq) + _____Na₃PO₄(aq)
$$\rightarrow$$
 _____Ca₃(PO₄)₂(s) + _____NaCl(aq)
B. _____BaCl₂(aq) + _____Na₂SO₄(aq) \rightarrow _____BaSO₄(s) + _____NaCl(aq)
C. _____FeCl₃(aq) + _____KSCN(aq) \rightarrow _____Fe(SCN)₃(aq) + _____KCl(aq)

Reaction type: _____

V. Sodium Carbonate and Hydrochloric Acid

	ons:				
Initial:	NaCO ₃	I	HCl		
Reaction	on with HCl:				
Splint:					
Balance: _	$\underline{\qquad}$ Na ₂ CO ₃ (s) + $\underline{\qquad}$	$_$ HCl(aq) \rightarrow $_$	_NaCl(aq) +	CO ₂ (g) +	$_H_2O($
* The read	ction can be considered	l as o ccurri ng in t	wo steps:		
Step 1:	$Na_2CO_3(s) + HCl(aq)$	→ NaC	$I(aq) + H_2CO_3(a)$	q)	
Step 2:	$H_2CO_3(aq)$	• $CO_2(g) + H_2O(g)$	1)		
Reaction t	ype: Step 1:		Step 2:		
OST-LAB (QUESTIONS				
1. What	evidence that a chemic	cal reaction is takin	ng place might ye	ou observe when:	
a.	Bananas ripen				
b.	Milk turns sour			_	
C.	An Alka-Seltzer table	et is added to wate	r		
d.	Burning a candle				
	te the following equations the following equations the second secon			decomposition, s	ingle

a. ____Al(s) + ____HBr(aq) \rightarrow ____AlBr₃(s) + ____H2(g) _____ b. ____H2(g) + ___Cl_2(g) \rightarrow ____HCl(g) _____ c. ___P4(s) + ___O2(g) \rightarrow ___P4O₁₀(s) _____

- 3. Write a <u>balanced</u> chemical equation for the following reactions write the formulas of the reactants and products, then balance. Identify the reaction type.
 - a. Ammonium nitrate forms nitrogen, oxygen and water.

Equation: _____ Reaction Type: _____

- b. Zinc chloride and sodium hydroxide form sodium chloride and zinc hydroxide.

Reaction Type: _____

4. Predict the products, then **balance** the chemical equations. Identify the type of reaction.

Equation:

Reaction Type

a. _Al + _O₂ \rightarrow _____ b. _AgCl + _Zn \rightarrow _____ + ___ c. _Pb(NO₃)₂ + _KI \rightarrow _____ + ___